



ARTICLE



Improvements in the Physical Properties of Agglomerated Milk Protein Isolate/Skim Milk Powder Mixtures Via Fluidized Bed Agglomeration

Chan Won Seo*

R&D Center, Seoul Dairy Cooperative, Ansan, Korea

Received: September 4, 2022
Revised: September 26, 2022
Accepted: September 26, 2022

*Corresponding author :
Chan Won Seo
R&D Center, Seoul Dairy Cooperative,
Ansan, Korea
Tel : +82-31-481-0146
Fax : +82-31-491-9179
E-mail : cwseo87@seoulmilk.co.kr

Copyright © 2022 Korean Society of Dairy Science and Biotechnology. This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<http://creativecommons.org/licenses/by-nc/3.0>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

ORCID

Chan Won Seo
<https://orcid.org/0000-0003-0787-2007>

Abstract

Protein-enriched dairy powder is widely consumed to promote muscle synthesis. Recently, in Korea, elderly people have also begun consuming protein powder products to prevent muscle loss. However, these protein-enriched powders have poor flowability and hydration properties because of the fine particles of spray-dried milk protein powder. Therefore, in this study, the fluidized bed agglomeration process was used to solve these problems. The rheological and physical properties of milk protein isolate (MPI)/skim milk powder (SMP) mixtures were effectively improved via fluidized bed agglomeration. The particle size of the MPI/SMP mixtures significantly increased from 35.7–58 μm to 118–136 μm , the flowability level improved from fair (21.4–26.3) to good (15.7–16.3), and the cohesiveness level changed from intermediate (1.27–1.36) to low (1.18–1.19) after fluidized bed agglomeration. In addition, the wetting time of the agglomerated MPI/SMP mixtures was effectively reduced to 4.67–58.3 s by fluidized bed agglomeration. These findings may be useful for manufacturing protein-enriched dairy powders with good instant properties.

Keywords

milk protein isolate, dairy powder, fluidized bed agglomeration, flowability, wettability

Introduction

High-protein foods are widely consumed for health purposes, and milk protein is an ideal ingredient owing to its nutritional benefits. Milk protein isolate (MPI) is a representative milk protein powder with a protein level above 90% on a dry matter basis that can provide all essential amino acids and valuable minerals, such as calcium and magnesium [1,2]. MPI contains both slow-acting protein (casein) and fast-acting protein (whey), which can supply sustained amino acids. Owing to these advantages, the demand for MPI is rapidly increasing with the growth of the high-protein food market [2,3]. MPI is generally manufactured by spray drying after removing lactose and mineral salts via ultrafiltration [1,4]. However, spray-dried MPI powder has poor flowability and wettability owing to its fine particle size and dense structure, which can lead to slow rehydration and lump formation [5,6].

Fluidized bed agglomeration can be defined as the size enlargement process from small fine particles to agglomerated large particles that is mainly used to improve the physical properties, such as flowability and wettability, of spray-dried powders with fine particles [7,8]. Currently, high-protein dairy powders are consumed in Korea to promote muscle synthesis and prevent muscle loss, and are generally manufactured by blending spray-dried protein and dairy-based powders with some nutritional ingredients, such as vitamins and minerals [9,10]. Thus, fluidized bed agglomeration is required to improve

the physical properties of high-protein dairy powders produced using the dry blending method. Some researchers have studied the effects of the fluidized bed agglomeration on dairy powders such as MPI and skim milk powder (SMP). Ji et al. [5] reported the improved wettability of MPI powders agglomerated by fluidized bed agglomeration, and Turchiuli et al. [11] observed the agglomerate growth of SMP particles during the fluid bed agglomeration process. However, the effect of the fluidized bed agglomeration process in the mixed system of MPI and SMP has not yet been studied. Therefore, the aim of this study was to investigate the effects of fluidized bed agglomeration on mixtures of MPI and SMP at different MPI/SMP ratios.

Materials and Methods

1. Materials

MPI (1900; 85% protein) was supplied by Vitalus Nutrition (Abbotsford, Canada), and SMP containing 37% protein and 49% lactose was obtained from Seoul Dairy Cooperative (Ansan, Korea). The mixtures of MPI and SMP were blended at ratios of 0:10 (SMP), 1:9 (MPI1/SMP9), 2:8 (MPI2/SMP8), and 3:7 (MPI3/SMP7) using a tumbler mixer (Bio-Components Inversina 2 L, Switzerland) for 30 min.

2. Fluidized bed agglomeration process

MPI/SMP mixture (1 kg) was placed in a product chamber and agglomeration was performed using a laboratory-scale fluidized bed granulator (Fluid Bed Lap System; Einsystem, Korea). Distilled water (binder solution) was pumped into the spray nozzle at a flow rate of 10 mL/min and sprayed at a pressure of 1.5 bar, and the product temperature was maintained at $55 \pm 1^\circ\text{C}$ during agglomeration. After spraying with distilled water for 40 min, the agglomerated powders were dried in hot air at 55°C for 10 min.

3. Measurement of particle size distribution

The particle size distributions of the fine and agglomerated MPI/SMP mixtures were determined using a laser diffraction particle size analyzer (Mastersizer 3000; Malvern Instruments, UK). The $D[3,2]$ (surface-weighted mean) and span (dispersion index) values were used to interpret the particle size distribution of the agglomerated MPI/SMP mixtures, calculated using Equations (1) and (2), respectively:

$$D[3,2] = \frac{\sum n_i d_i^3}{\sum n_i d_i^2} \quad (1)$$

$$\text{Span} = \frac{Dv90 - Dv10}{Dv50} \quad (2)$$

where d_i is the average diameter of the i -th class interval and n_i is the number of

particles with diameter d_i . $Dv10$, $Dv50$, and $Dv90$ are the average particle sizes at 10%, 50%, and 90% of the cumulative size distribution, respectively.

4. Density measurement

Bulk density (ρ_{bulk}) and tapped density (ρ_{tapped}) of the MPI/SMP mixtures were measured by the method of Wu et al. [6] with slight modification using a tap density meter (BT-301; Battersize Instrument, China). The ρ_{bulk} and ρ_{tapped} values were calculated as the mass/volume ratios of the powder samples before and after 1250 taps, respectively.

5. Flowability and cohesiveness measurements

The flowability and cohesiveness of the MPI/SMP mixtures were evaluated using the Carr index (CI) and the Hausner ratio (HR) [7]. The CI and HR values were calculated using ρ_{bulk} and ρ_{tapped} using Equations (3) and (4), respectively:

$$CI = \frac{\rho_{tapped} - \rho_{bulk}}{\rho_{tapped}} \times 100 \quad (3)$$

$$HR = \frac{\rho_{tapped}}{\rho_{bulk}} \quad (4)$$

The classification of powder flowability (%) and cohesiveness based on CI and HR values are presented in Table 1.

6. Measurement of wetting time

The wetting time (t_w) was measured to evaluate the rehydration properties of fine and agglomerated MPI/SMP mixtures. According to the method of Atalar and Yazici [7], the wetting time was determined as the time required to complete wetting and immersion of 5 g of powder in 100 mL of distilled water ($24 \pm 1^\circ\text{C}$).

Table 1. Classification of powder flowability and cohesiveness based on Carr index (CI) and Hausner ratio (HR)

		Classification
CI (flowability)	<15	Very good
	15–20	Good
	20–35	Fair
	35–45	Bad
	>45	Very bad
HR (cohesiveness)	<1.2	Low
	1.2–1.4	Intermediate
	>1.4	High

CI, Carr index; HR, Hausner ratio.

7. Statistical analysis

Statistical analysis was conducted using one-way analysis of variance, followed by Tukey's test, using the IBM SPSS Statistics 24 software suite (IBM Software, USA). The level of statistical significance was set at $p < 0.05$.

Results

1. Particle size distribution

The particle size distribution profiles of the fine and agglomerated MPI/SMP mixtures with different MPI ratios are shown in Fig. 1. The fine powder samples exhibited monomodal distributions, which moved to a large region after the fluidized bed agglomeration process. These results indicate that the particle size of the powders increased as they agglomerated during the fluidized bed agglomeration process.

The average particle diameter and distribution of the powder were characterized using the $D[3,2]$ (surface-weighted mean) and span values, respectively. As shown in Table 2, the $D[3,2]$ value of fine SMP was 58.0 μm , which decreased to 35.7 μm (MPI3/SMP7) with an increase in the MPI ratio due to the small particle size of fine MPI (22.7 μm). After the fluidized bed agglomeration process, the $D[3,2]$ values of the MPI/SMP mixtures increased. The span values of the fine MPI/SMP mixtures increased from 1.48 to 1.90 as the MPI ratio increased due to the high span value of fine MPI




Fig. 1. Particle size distribution of fine and agglomerated MPI/SMP mixtures with different MPI ratios: (A) SMP (0:10), (B) MPI1/SMP9 (1:9), (C) MPI2/SMP8 (2:8), and (D) MPI3/SMP7 (3:7). MPI, milk protein isolate; SMP, skim milk powder.

Table 2. Particle size and distribution of fine and agglomerated MPI/SMP mixtures with different MPI ratios

	Sample	D[3,2] (μm)	Span
Fine powder	MPI	22.7 \pm 0.40 ^a	2.07 \pm 0.02 ^a
	SMP	58.0 \pm 0.85 ^b	1.48 \pm 0.01 ^b
	MPI1/SMP9	47.8 \pm 1.26 ^c	1.67 \pm 0.01 ^c
	MPI2/SMP8	40.4 \pm 0.15 ^d	1.79 \pm 0.01 ^{df}
	MPI3/SMP7	35.7 \pm 0.15 ^d	1.90 \pm 0.01 ^e
Agglomerated powder	SMP	151 \pm 3.79 ^e	1.81 \pm 0.03 ^{df}
	MPI1/SMP9	136 \pm 1.00 ^f	1.82 \pm 0.04 ^f
	MPI2/SMP8	127 \pm 1.73 ^g	1.75 \pm 0.02 ^{df}
	MPI3/SMP7	118 \pm 1.15 ^h	1.75 \pm 0.03 ^d

Values are the means of three measurements \pm SD.

MPI (milk protein isolate)/SMP (skim milk powder) mixtures were mixed at different ratios of 1:9 (MPI1/SMP9), 2:8 (MPI2/SMP8), and 3:7 (MPI3/SMP7).

^{a-h} Mean values in the same line with different letters are significantly different ($p < 0.05$).

(2.07), indicating a broad size distribution and low homogeneity. Except for the agglomerated MPI1/SMP9 mixture, the span values of the agglomerated MPI/SMP mixtures were lower than those of the non-agglomerated MPI/SMP mixtures, indicating a narrow size distribution with homogenous powders.

2. Bulk and tapped density

The densities of the fine and agglomerated MPI/SMP mixtures with different MPI ratios are shown in Fig. 2. The bulk density (ρ_{bulk}) and tapped density (ρ_{tapped}) values of the fine powder samples were in the range 0.40–0.47 and 0.54–0.58 g/mL, respectively. After the fluidized bed agglomeration process, ρ_{bulk} and ρ_{tapped} values decreased to 0.26–0.32 and 0.31–0.38 g/mL, respectively, and decreased with an increase in the MPI ratio.

3. Flowability and cohesiveness

The flowability and cohesiveness of food powders can be characterized using CI and HR values, respectively. As shown in Table 3, the CI value of the fine SMP sample showed good flowability (CI: 15–20) and increased from 18.4 to 26.3 as the MPI ratio increased, indicating fair flowability (CI: 20–35). The HR value of the fine SMP sample was 1.23 and that of the fine MPI/SMP mixtures increased with increasing MPI ratio, indicating intermediate cohesiveness (HR: 1.2–1.4). In contrast, the CI and HR values of the agglomerated MPI/SMP mixtures showed good flowability (CI<15) and low cohesiveness (HR<1.2). These results indicate that the physical properties of MPI/SMP mixtures can be effectively improved via fluidized bed agglomeration.

4. Wetting time

The wettabilities of the fine and agglomerated MPI/SMP mixtures with different MPI ratios are listed in Table 3. The fine powder samples exhibited a long wetting time (t_w) >120 s and floated on the water surface with lump formation. In contrast, the wetting




Fig. 2. Densities of fine and agglomerated MPI/SMP mixtures with different MPI ratios: (A) bulk density (ρ_{bulk}), (B) tapped density (ρ_{tapped}). MPI, milk protein isolate; SMP, skim milk powder.

Table 3. Physical and hydration properties of fine and agglomerated MPI/SMP mixtures with different MPI ratios

	Sample	CI (%)	HR	t_w (s)
Fine powder	SMP	18.4±1.05 ^a	1.23±0.02 ^a	>120
	MPI1/SMP9	21.4±1.06 ^b	1.27±0.02 ^b	>120
	MPI2/SMP8	24.3±1.60 ^c	1.32±0.03 ^c	>120
	MPI3/SMP7	26.3±0.49 ^c	1.36±0.01 ^c	>120
Agglomerated powder	SMP	18.4±0.22 ^a	1.23±0.01 ^a	4.67±1.15 ^a
	MPI1/SMP9	16.3±0.46 ^{ad}	1.19±0.01 ^a	10.3±1.52 ^b
	MPI2/SMP8	15.5±1.28 ^d	1.18±0.02 ^a	16.7±3.06 ^c
	MPI3/SMP7	15.7±0.89 ^{ad}	1.19±0.01 ^a	58.3±11.7 ^d

Values are the means of three measurements±SD.

MPI (milk protein isolate)/SMP (skim milk powder) mixtures were mixed at different ratios of 1:9 (MPI1/SMP9), 2:8 (MPI2/SMP8), and 3:7 (MPI3/SMP7).

^{a-d} Mean values in the same line with different letters are significantly different ($p<0.05$).

CI, Carr index; HR, Hausner ratio.

time of agglomerated powder samples significantly decreased to 4.67–58.3 s. With the exception of the MPI3/SMP7 sample, agglomerated powder samples exhibited a faster wetting time than the instant criterion (<20 s) described by Chever et al. [12].

Discussion

Particle size distribution is an important property that is closely related to other physical properties, such as density, flowability, and wettability, of powdered products. As evidenced by the increased particle size, the fine powder samples were effectively agglomerated by the fluidized bed process. This agglomeration of fine powders can be attributed to the formation of solid bridges between fine particles during the fluidized bed process. In general, the fluidized bed agglomeration process is divided into three steps: wetting, coalescence, and consolidation. The fine powders were fluidized by hot air and their surfaces were wet by the sprayed binder solution. The wet powder particles collided with each other and bound together, leading to the formation of liquid bridges. Subsequently, the liquid bridges become solid bridges by hot air, which leads to the consolidation of agglomerated particles [13,14]. Therefore, the fluidized bed process can effectively increase the particle size and produce agglomerated powders.

The main purpose of the fluidized bed agglomeration process is to achieve uniform particle size and distribution. Span value is generally used to evaluate the homogeneity and polydispersity of food powders. The high span value of the agglomerated SMP sample may be attributed to the high friability of the agglomerated SMP particles owing to friction within the fluidized bed [15,16]. In contrast, the span value of the agglomerated MPI3/SMP7 sample was significantly lower than that of the non-agglomerated MPI3/SMP7 sample. This reduction in the span value indicates that the broad distribution of the fine powders was transformed into a narrow size distribution with homogenous powders via fluidized bed agglomeration. From these results, it can be concluded that the particle size distribution of the MPI/SMP mixtures can be effectively improved using the fluidized bed agglomeration process.

Fluidized bed agglomeration led to a lower density than that of the non-agglomerated powder samples. This result can be attributed to the increased particle size and porous structure of the agglomerated powders, which can cause additional voids and result in a lower density [5,12]. In contrast, the fine powder samples showed higher densities owing to their dense structures. Consequently, the decreased ρ_{bulk} and ρ_{tapped} values of the agglomerated MPI/SMP mixtures can be explained by the additional voids within the porous structure of the agglomerated powders.

CI and HR values of the MPI/SMP mixtures were significantly improved via fluidized bed agglomeration. The improved flowability and cohesiveness of the agglomerated powders can be explained by the increased particle size, which can decrease the friction between the particles owing to the reduction in the particle surface area [14,17]. Fluidized bed agglomeration also improved the wettability of the MPI/SMP mixtures. This improved wettability can be attributed to the porous structure of the agglomerated powder. The large porous structure allows water to easily penetrate into the void spaces of the agglomerated particles, facilitating the wetting of the powders [6,8]. These results suggest that the physical properties, such as flowability and wettability, of the MPI/SMP mixtures can be effectively improved to obtain increased particle size and porous structure via the fluidized bed agglomeration process.

Conclusion

The high-protein dairy powder is being consumed to promote muscle synthesis and prevent muscle loss. However, spray-dried milk proteins have poor flowability and wettability owing to their fine particle sizes and high inter-particle cohesive forces. To resolve these problems, this study investigated the effect of fluidized bed agglomeration on MPI/SMP mixtures with varying MPI/SMP ratios. The fluidized bed process increased the particle size from 35.7–58 μm to 118–136 μm , and both CI and HR values of MPI/SMP mixtures decreased from 21.4–26.3 and 1.27–1.36 to 15.7–16.3 and 1.18–1.19, respectively. These results showed that the flowability level was improved from fair (20–35) to good (15–20) and cohesiveness level was also changed from intermediate (1.2–1.4) to low (<1.2) after the fluidized bed agglomeration. Moreover, the wetting time of agglomerated MPI/SMP mixtures was effectively reduced to 4.67–58.3 s by the fluidized bed agglomeration. In particular, agglomerated MPI1/SMP9 and MPI2/SMP8 samples exhibited a faster wetting time than the instant criterion (<20 s). The results of this study may be useful in manufacturing protein-enriched dairy powders with good instant properties.

Conflict of Interest

The authors declare no potential conflict of interest.

References

1. Uluko H, Liu L, Lv JP, Zhang SW. Functional characteristics of milk protein concentrates and their modification. *Crit Rev Food Sci Nutr*. 2016;56:1193–1208.
2. Meena GS, Singh AK, Panjagari NR, Arora S. Milk protein concentrates: opportunities and challenges. *J Food Sci Technol*. 2017;54:3010–3024.
3. Agarwal S, Beausire RLW, Patel S, Patel H. Innovative uses of milk protein concentrates in product development. *J Food Sci*. 2015;80:A23–A29.
4. Corredig M, Nair PK, Li Y, Eshpari H, Zhao Z. Invited review: understanding the behavior of caseins in milk concentrates. *J Dairy Sci*. 2019;102:4772–4782.
5. Ji J, Cronin K, Fitzpatrick J, Fenelon M, Miao S. Effects of fluid bed agglomeration on the structure modification and reconstitution behaviour of milk protein isolate powders. *J Food Eng*. 2015;167:175–182.
6. Wu S, Fitzpatrick J, Cronin K, Maidannyk V, Miao S. Effects of spraying surfactants in a fluidised bed on the rehydration behaviour of milk protein isolate powder. *J Food Eng*. 2020;109694.
7. Atalar I, Yazici F. Influence of top spray fluidized bed agglomeration conditions on the reconstitution property and structure modification of skim yoghurt powder. *J Food Process Preserv*. 2017;42:e13414.
8. Custodio GR, de Souza LFG, Nitz M, Andreola K. A protein powder agglomeration process using açai pulp as the binder: an analysis of the process parameters. *Adv*

- Powder Technol. 2020;31:3551-3561.
9. Jang JH, Oh N. Current trends in the food industry for development of milk protein based products. *Livest Food Sci Ind.* 2021;10:75-84.
 10. Kang Y, Kim N, Choi YJ, Lee Y, Yun J, Park SJ, et al. Leucine-enriched protein supplementation increases lean body mass in healthy Korean adults aged 50 years and older: a randomized, double-blind, placebo controlled trial. *Nutrients.* 2020;12:1816.
 11. Turchiuli C, Smail R, Dumoulin E. Fluidized bed agglomeration of skim milk powder: analysis of sampling for the follow-up of agglomerate growth. *Powder Technol.* 2013;238:161-168.
 12. Chever S, Méjean S, Dolivet A, Mei F, Den Boer CM, Le Barzic G, et al. Agglomeration during spray drying: physical and rehydration properties of whole milk/sugar mixture powders. *LWT-Food Sci Technol.* 2017;83:33-41.
 13. Atalar I, Yazici F. Effect of different binders on reconstitution behaviors and physical, structural, and morphological properties of fluidized bed agglomerated yoghurt powder. *Dry Technol.* 2019;37:1656-1664.
 14. Lim DH, Letona A, Lee M, Lim D, Han NS, Chung D. Fluidized-bed granulation of probiotics-encapsulated spray-dried skim milk powder: effects of a fluidizing aid, moisture-activation and dehydration. *Foods.* 2021;10:1600.
 15. Lee D, Min G, Roh W, Yoo B. Effect of various types of sugar binder on the physical properties of gum powders prepared via fluidized-bed agglomeration. *Foods.* 2021;10:1387.
 16. Nascimento RF, Rosa JG, Andreola K, Taranto OP. Wettability improvement of pea protein isolate agglomerated in pulsed fluid bed. *Part Sci Technol.* 2020;38:511-521.
 17. Park J, Yoo B. Particle agglomeration of gum mixture thickeners used for dysphagia diets. *J Food Eng.* 2020;279:109958.